


PROTECT[®] AND CARE



Colour: black

PRODUCT FILE

Nitrile gloves - disposable - powder free Black 5g

Ref. A00137

Specifications

- Single use
- Suitable for HACCP and medical purposes
- Prevent contamination & protect against viruses, bacteria & fungi
- 100% Nitrile: stronger than other glove materials such as latex, natural rubber or vinyl.
- Are flexible and offer much greater tactile sensitivity than standard diagnostic gloves.
- Good puncture and tear resistance
- Hypo allergenic: - Latex free

- Powder free (less chance of irritation.)


- Impermeable when used in wet or greasy environments
- Comfortable fit:
 - fitted with roll cuffs for easy donning and doffing
 - ambidexterous
- The roughened fingertips ensure a secure grip, even in wet environments

- Sizes:

	Length mm	Palm width mm
S 6-7	240	80 +/-10
M 7-8	240	95 +/-10
L 8-9	240	110 +/-10
XL 9-10	240	+110 +/-10



- Clear size indication with color
- Packaging with barcode and QR code
- Info in 27 languages: NL, FR, EN, DE, ES, IT, PT, BG, CS, DA, EL, FI, HU, HR, KK, LT, LV, NO, PL, RO, RU, SK, SL, SR BS, SV, TR & UK



Packaging

Dispenser box of 100 pcs/box - 10 boxes/master box = 1000 gloves (500 pairs)

EN ISO 374-1:2016
+A1:2018
TYPE B

GESCHIKT VOOR HACCP EN MEDISCHE DOELEINDEN
APPROPRIE POUR HACCP ET SOINS MEDICALES
SUITABLE FOR HACCP AND MEDICAL PURPOSES
FÜR HACCP UND MEDIZINISCHE ZWECKE GEEIGNET

Letter	Chemical	CAS	Permeation Level EN ISO 374-1:2016+A1:2018 Type B	Degradation Level EN ISO 374-4:2019
K	40% Sodium Hydroxide	1310-73-2	6	-8.3%
P	30% Hydrogen Peroxide	7722-84-1	2	34.1%
T	37% Formaldehyde	50-00-0	4	34.3%

Permeation Performance Level

Measured breakthrough time (minutes)	1	2	3	4	5	6
	>10	>30	>60	>120	>240	>480

In compliance with: MDR (EU) 2017/745 Class I
PPE Regulation (EU) 2016/425: Category III
EN ISO 21420:2020
EN 455-1-2020, EN 455-2-2015, EN 455-3-2015
ASTM D 6978, ASTM F 1671

EU Type Examination and Module C2 Notified Body 2777
SATRA TECHNOLOGY Europe Limited, Bracetown Business Park,
Clonee Dublin 15, D15 YN2P Republic of Ireland

Test acc. to ASTM D 6978 - Test Chemotherapy

TEST CHEMOTHERAPY DRUGS AND CONCENTRATION	MINIMUM BREAKTHROUGH DETECTION TIME (Specimen 1/2/3) (Minutes)
Carmustine (BCNU) 3.3 mg/ml (3,300 ppm)	14.7 (15.1, 14.2, 15.8)
Cisplatin 1.0 mg/ml (1,000 ppm)	240 min.
Cyclophosphamide (Cytozen) 20 mg/ml (20,000 ppm)	240 min.
Cytarabine 100 mg/ml (100,000 ppm)	240 min.
Doxorubicine (DUC) 10.0 mg/ml (10,000 ppm)	240 min.
Doxorubicin Hydrochloride 2.0 mg/ml (2,000 ppm)	240 min.
Etoposide (Etopos) 20.0 mg/ml (20,000 ppm)	240 min.
Fluorouracil 50.0 mg/ml (50,000 ppm)	240 min.
Irinotecan 50.0 mg/ml (50,000 ppm)	240 min.
Methotrexate 25 mg/ml (25,000 ppm)	240 min.
Mitomycin C 0.5 mg/ml (500 ppm)	240 min.
Misocarbonyl 2.0 mg/ml (2,000 ppm)	240 min.
Paclitaxel (Taxol) 5.0 mg/ml (5,000 ppm)	240 min.
Procarbazine 10.0 mg/ml (10,000 ppm)	58.8 (110.0, 58.8, 67.0)
Vincristine Sulfate 1.0 mg/ml (1,000 ppm)	240 min.

ISO 374-5:2016

TESTED IN ACCORDANCE TO:
EN ISO 374-5:2016 and EN ISO 374-4:2019
Degradation results indicate the change in puncture resistance of the gloves after exposure to the challenge chemical
PROTECTION AGAINST BACTERIA, VIRUSES & FUNGI- PASS