


PROTECT[®] AND CARE



PRODUCT FILE

Nitrile medical examination gloves 5 gsm

Non sterile

Ref. A00204

Specifications

- Single use
- Suitable for HACCP and medical purposes
- Medical and protective: prevent contamination and protect against viruses, bacteria & fungi
- 100% Nitrile: stronger than other glove materials such as latex, natural rubber or vinyl.
- Are flexible and offer much greater tactile sensitivity than standard diagnostic gloves.
- Good puncture and tear resistance
- Hypo allergenic: Latex and powder- free
- Impervious when used in wet or greasy working environment
- Comfortable fit:
 - fitted with roll cuffs for easy donning and doffing
 - ambidextrous
- The roughened fingertips ensure a secure grip, even in wet environments

Sizes:

	Length mm	Palm width mm
S 6-7	240	80 +/-10
M 7-8	240	95 +/-10
L 8-9	240	105 +/-10
XL 9-10	240	115 +/-5

- Clear size indication with color
- Packaging with barcode and QR code
- Info in 27 languages: NL, FR, EN, DE, ES, IT, PT, BG, CS, DA, EL, FI, HU, HR, KK, LT, LV, NO, PL, RO, RU, SK, SL, SR, BS, SV, TR & UK



Packaging

Dispenser box of 100 pcs/box 10 boxes/master box = 1000 gloves (500 pairs)

Certificatie

EN ISO 374-1:2016
+A1:2018
TYPE B



Letter	Chemical	CAS	Permeation Level EN ISO 374-1:2016+A1:2018 Type B	Degradation Level EN ISO 374-4:2019
K	40% Sodium Hydroxide	1310-73-2	6	-8.3%
P	30% Hydrogen Peroxide	7722-84-1	2	34.1%
T	37% Formaldehyde	50-00-0	4	34.3%

Permeation Performance Level	1	2	3	4	5	6
Measured breakthrough time (minutes)	>10	>30	>60	>120	>240	>480

In compliance with: MDR (EU) 2017/745 Class I
PPE Regulation (EU) 2016/425: Category III
EN ISO 21420:2020
EN 455-1-2020, EN 455-2-2015, EN 455-3-2015
EU Type Examination and Module C2 Notified Body 2777
SATRA TECHNOLOGY Europe Limited, Bracetown Business Park,
Clonlee Dublin 15, D15 YN2P, Republic of Ireland

ISO 374-5:2016



TESTED IN ACCORDANCE TO:
EN ISO 374-5:2016 and EN ISO 374-4:2019
Degradation results indicate the change in puncture resistance
of the gloves after exposure to the challenge chemical
PROTECTION AGAINST BACTERIA, VIRUSES & FUNGI- PASS

GESCHIKT VOOR HACCP EN MEDISCHE DOELEINDEN
APPROPRIE POUR HACCP ET SOINS MEDICALES
SUITABLE FOR HACCP AND MEDICAL PURPOSES
FÜR HACCP UND MEDIZINISCHE ZWECKE GEEIGNET



SPECIFICATION OF NITRILE DISPOSABLE GLOVES powder-free						
Size Spec		Small	Medium	Large	X-large	Tolerance
Total Length (mm)		240	240	240	240	min
Width of palm (mm)		85	95	105	115	± 5
Thickness (mm)	Finger	0.09	0.09	0.09	0.09	± 0.02
	Palm	0. 07	0.07	0.07	0.07	± 0.02
	Wrist	0.05	0.05	0.05	0.05	± 0.02
Weight per piece (gram)		3.0	3.5	3.9	4.3	± 0.3
T ensile-strength (≧)min		14Mpa	14Mpa	14Mpa	14MPa	min
Elongation (%)min	Before aging:	500	500	500	500	min
	After aging:	400	400	400	400	min